Velocity^{ALPHA} Highspeed Surgical Drill System

Velocity^{ALPHA}

HIGHSPEED SURGICAL DRILL SYSTEM

*30% less weight**

*25% more power**

*25% slimmer**

Minimal thermal warming

**compared to previous Velocity generation*

Based on more than 30 years of highspeed motordrill experience and continuous development, adeor[®] launches the all new electric Velocity[®] Alpha highspeed surgical drill system.

The Velocity[®] Alpha system offers distinct advantages to all parties: Reliable performance and unique surgical features for the surgeon, easiest setup and handling for OR staff, significantly reduced cost and cash-flow requirements for purchasing departments.

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VELOCITY

The universal cutter system

80% less inventory cost / different cutter models


No more cutter identification hassles

Adjustable cutter exposure

Easiest handling

UNIVERSAL
CUTTER

Velocity[®] Alpha universal cutter system

> 80% less cutter codes and inventory cost

Instead of more than 200 to 300 different cutters and product codes, Velocity Alpha only has 50 codes without compromising variety. This means less stock requirement and free cash-flow for the hospital.

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> no more cutter identification

Warehouse managers, nurses, surgeons, purchasing officers will all benefit from the simplicity of the universal cutter design. Less codes mean less hassle.

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
> Improved safety

The common risk of using the wrong cutter with the wrong nosepiece is eliminated.

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
UNIVERSAL
CUTTER

One universal cutter for all straight and angled Velocity Alpha nosepieces.

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One cutter - 8 exposure positions.
Peerless flexibility.

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
> variable cutter exposure

Depending on the procedure, surgical area and other factors, surgeons want cutters to protrude more or less from the nosepiece. Now, they can select from 8 different positions for each tool (exposure range = 8 mm).

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ONE-CLICK
COUPLING

One click coupling mechanism.
No twisting, no pulling, no headache.

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> one-click coupling

Inserting, removing or changing cutters has never been easier. To insert the tool, simply push the cutter into the nosepiece until you hear and feel the "click". That's it. Less training, less risk of errors.

OR Time is expensive. The automatic one-click coupling system saves significant time, during set-up as well as intra-operative.

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One click to connect any cutter

One click to connect any attachment

Straight nosepieces

> Super slim - packed with unique features.

Velocity Alpha straight nosepieces are designed with the toughest challenges in microsurgical procedures in mind. They are built to last and offer reliable performance without excessive heat, vibration or noise. All Velocity Alpha nosepieces offer a unique set of features:

- > Universal cutter system
- > Variable cutter exposure
- > One click coupling
- > Super slim shaft diameter

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total length 170 mm

Angled nosepieces

> Super slim - packed with unique features.

Velocity Alpha angled nosepieces are designed with the toughest challenges in microsurgical procedures in mind. They are built to last and offer reliable performance without excessive heat, vibration or noise. All Velocity Alpha nosepieces offer a unique set of features:

- > Universal cutter system
- > Variable cutter exposure
- > One click coupling
- > Super slim shaft diameter
- > 25° angulation for enhanced vision

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VMN-A1

short attachment angled
working length 40 mm
total length 120 mm

VMN-A2

medium attachment angled
working length 70 mm
total length 150 mm

VMN-A3

long attachment angled
working length 100 mm
total length 180 mm

VMN-A4

xl attachment angled
working length 130 mm
total length 210 mm

VMN-A5

xxl attachment angled
working length 170 mm
total length 250 mm

VMN-A6

xxxl attachment angled
working length 220 mm
total length 300 mm

Craniotomy attachments

> Modular setup

The Velocity Alpha craniotomy attachments consist of a combination of a solid heavy duty craniotome base which can be coupled with six different duraguards. The duraguards are available in three lengths, fixed or rotatable. The duraguard design facilitates smooth and fast cutting:

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VMN-C1

craniotome base

holds fixed and rotatable duraguards
Safety lock mechanism

VMCH-1

short duraguard

12 mm cutting width
fixed version

VMCH-2

medium duraguard

15 mm cutting width
fixed version

VMCH-3

long duraguard

22 mm cutting width
fixed version

VMCH-1R

short duraguard

12 mm cutting width
rotatable version

VMCH-2R

medium duraguard

15 mm cutting width
rotatable version

VMCH-3R

long duraguard

22 mm cutting width
rotatable version

Speed reducer

> Ergonomic shape. High torque.

The Velocity Alpha speed reducer is designed for fast and efficient cranial perforation. For best performance please use adeor cranial perforators (see pages 20-21).

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Highspeed motor

> High speed and torque - minimal thermal warming.

The Velocity highspeed motor is an ideal combination of ergonomic shape, power and inbuilt insulation. Even prolonged surgical procedures are possible without the need for intra-operative cooling periods.

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VMT-102E

Speed reducer

attaches to all cranial perforators with hudson fitting
high torque / 1:64 ratio

VMM-100EC

80.000 rpm highspeed motor

with integrated 3.5 m cable
autoclavable < 134°C

Velocity® Control Unit

The compact control unit has been designed utilising feedback from numerous surgeons in various regions of the world. The result is a console that is easy to operate with only a minimum of buttons. No lengthy menus to navigate, no overload of functionality. This reduces the likelihood of operating errors, time needed for installation and set-up as well as training requirements.

The control unit meets or exceeds all applicable safety standards. Its dimensions are a new benchmark for high-speed drill systems. It can be easily placed close to the user and almost without any hindrance to other medical apparatus within the OR.

Most functions of the control unit can also be operated by the surgeon using the foot switch. This is not just time-saving but also reduces the need to communicate between sterile and non-sterile OR personnel.

The control unit also features an inbuilt irrigation pump with adjustable flow rate. This feature eliminates the need for an additional irrigation system for the drilling process.

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**VMM-EC**

Control unit with integrated irrigation pump

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FOOT CONTROL

- > Large, easy to find foot print
- > Wide-angle accessibility
- > 3 colour-coded buttons to control rotation, operation mode and irrigation
- > Accelerator or On/Off Mode
- > Safety guard

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
Custom stainless steel sterilization tray

- > secure reprocessing and storage
- > customized inlays for proper fit of all components
- > pictograms for ease of use

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Velocity[®] Alpha oil-spray

- > one universal spray for cleaning and lubrication
- > to be used with custom made adapter OA-100
- > extends life span of components

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PST-200

sterilization tray PPSU plastic
custom made inlays
autoclavable < 134°C

OA-100

oil-spray adapter
screw-on profile
autoclavable < 134°C

PST-201

sterilization tray stainless steel
custom made inlays
autoclavable < 134°C

OS-100

oil-spray
400 ml can
12 cans per box

Irrigation tubing set and adapters

- > customized irrigation adapters available for all nosepieces
- > can be connected to single-use sterile tubing set
- > targeted irrigation of drilling area
- > adjustable flow rate

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VMN-IR-1

short irrigation adapter

compatible with Velocity Alpha nosepieces
VMN-S1 and VMN-A1

VMN-IR-2

medium irrigation adapter

compatible with Velocity Alpha nosepieces
VMN-S2 and VMN-A2

VMN-IR-3

long irrigation adapter

compatible with Velocity Alpha nosepieces
VMN-S3 and VMN-A3

VMN-IR-4

xl irrigation adapter

compatible with Velocity Alpha nosepiece
VMN-A4

VMN-IR-5

xxl irrigation adapter

compatible with Velocity Alpha nosepiece
VMN-A5

VMN-IR-6

xxxl irrigation adapter

compatible with Velocity Alpha nosepiece
VMN-A6

Velocity® Alpha System overview

> Straight attachments

VMN-S1
40 mm

VMN-S2
70 mm

VMN-S3
100 mm

> Duraguards and craniotome base

VMCH-1
12 mm fixed

VMCH-1R
12 mm rotatable

VMCH-2
15 mm fixed

VMCH-2R
15 mm rotatable

VMCH-3
22 mm fixed

VMCH-3R
22 mm rotatable

VMN-C1
craniotome base

> Angled attachments

VMN-A1
40 mm

VMN-A3
100 mm

VMN-A5
170 mm

VMN-A2
70 mm

VMN-A4
130 mm

VMN-A6
220 mm

> Speed reducer

VMT-102E
1:64 ratio

> Irrigation tubing set & adapters

See page 15
for details

> Highspeed motor

VMM-100EC
80.000 rpm

> Oil adapter & Oil spray

OA-100 & OS-100
oil adapter & oil spray

> Control unit

VMM-EC-220
integrated irrigation

> Sterilization trays

PST-201
stainless steel sterilization tray

> Foot control

VMME-FS
water resistant

PST-200
PPSU sterilization tray

Velocity® Alpha universal cutter range

One universal cutter for all Velocity Alpha nosepieces.

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Meridian® Classic cranial perforator range

CRANIAL PERFORATORS

Unique Benefits

- › Popular design for maximized familiarity
- › Suitable for skull thicknesses as thin as 2 mm
- › Safe release mechanism
- › Suitable for even and uneven bone
- › Universal Hudson fitting
- › Durable cutting efficiency
- › Different housing colors to easily identify the different models
- › Available in three most commonly used sizes

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Single-Use Cranial Perforators
for use on skulls as thin as 2.0 mm

Meridian^{Classic} Single-Use Perforator
PER C6-9-S 6 / 9 mm (inner / outer diameter)
single-use, sterile
replaces Codman 26-1223

Meridian^{Classic} Single-Use Perforator
PER C8-11-S 8 / 11 mm (inner / outer diameter)
single-use, sterile
replaces Codman 26-1222

Meridian^{Classic} Single-Use Perforator
PER C11-14-S 11 / 14 mm (inner / outer diameter)
single-use, sterile
replaces Codman 26-1221

System Components / Systemkomponenten

Velocity Alpha
system components

| | |
|--|------------|
| Electric control unit Highspeed / Trepan 220V version | VMM-EC-220 |
| Electric control unit Highspeed / Trepan 110V version | VMM-EC-110 |
| Electric highspeed motor, 80.000 rpm, integrated cable 3.5 m | VMM-100EC |
| Electric foot switch | VMME-FS |
| Irrigation tubing set, 6 pcs, sterile | VMM-ITS-6 |
| Straight nosepiece short working length 40 mm | VMN-S1 |
| Straight nosepiece medium working length 70 mm | VMN-S2 |
| Straight nosepiece long working length 100 mm | VMN-S3 |
| Angled nosepiece short working length 40 mm | VMN-A1 |
| Angled nosepiece medium working length 70 mm | VMN-A2 |
| Angled nosepiece long working length 100 mm | VMN-A3 |
| Angled nosepiece XL working length 130 mm | VMN-A4 |
| Angled nosepiece XXL working length 170 mm | VMN-A5 |
| Angled nosepiece XXXL working length 220 mm | VMN-A6 |
| Irrigation attachment working length 40 mm | VMN-IR-1 |
| Irrigation attachment working length 70 mm | VMN-IR-2 |
| Irrigation attachment working length 100 mm | VMN-IR-3 |
| Irrigation attachment working length 130 mm | VMN-IR-4 |
| Irrigation attachment working length 170 mm | VMN-IR-5 |
| Irrigation attachment working length 220 mm | VMN-IR-6 |

System Components / Systemkomponenten

Velocity Alpha
system components

| | |
|--|----------|
| Craniotome base | VMN-C1 |
| Duraguard fixed, short 12 mm | VMCH-1 |
| Duraguard fixed, medium 15 mm | VMCH-2 |
| Duraguard fixed, long 22 mm | VMCH-3 |
| Duraguard rotatable, short 12 mm | VMCH-1R |
| Duraguard rotatable, medium 15 mm | VMCH-2R |
| Duraguard rotatable, long 22 mm | VMCH-3R |
| Trepan Hudson adapter, low speed, 1:64 ratio | VMT-102E |
| Highspeed oil spray, 12 x 400 ml | OS-100 |
| Custom made oil spray nozzle adapter, screw-on profile, reusable | OA-100 |
| Transport & Storage Case suitable for air transport | TSC-200 |
| Sterilization Tray plastic, customized | PST-200 |
| Sterilization Tray stainless steel, customized | PST-201 |

SPARE PARTS

| | |
|---------------------------------|-----------|
| Bar for footswitch | VMME-FS-H |
| Irrigation stand | VMM-EC-S |
| Irrigation tube clip set 5 pcs. | VMM-CLIPS |

Technical Data / Technische Daten

| | Velocity System EU Version 220 V | Velocity System US Version 110 V |
|--|---|---|
| Supply voltage | 220 - 240 V | 100 - 130 V |
| Permitted voltage fluctuation | ± 10 % | ± 10 % |
| Nominal current | 0.1 - 0.8 A | 0.2 - 1.7 A |
| Frequency | 50 - 60 Hz | 50 - 60 Hz |
| Mains fuse | 2 x 250 V - T 1.25 AH | 2 x 250 V - T 1.25 AH |
| Max. power consumption | 170 VA | 170 VA |
| Max. mechanical output power | 80 W | 80 W |
| Torque level | 2.0 ... 80 Ncm | 2.0 ... 80 Ncm |
| Motor speed range | 1.000 - 80.000 rpm | 1.000 - 80.000 rpm |
| Irrigation flow rate @ 65 / 80 / 100 % | 65 / 80 / 100 ml/min | 65 / 80 / 100 ml/min |
| Dimensions | 255 x 256 x 91 mm | 255 x 256 x 91 mm |
| Control unit weight | 2.3 kg | 2.3 kg |
| Protection class | II | II |
| Applied part | Type B | Type B |
| Integrated irrigation pump | yes | yes |
| Foot control / danger zone M / ANSI/AAMI | AP / IEC 60601-1 | AP / IEC 60601-1 |
| Water resistance of foot control / 1 m depth for 1 h | IPX8 / IEC 60529 | IPX8 / IEC 60529 |
| Pollution level | 2 | 2 |
| Overvoltage category | II | II |
| Max. altitude / above sea level | 3.000 m | 3.000 m |
| Short term storage temperature | 0°C ... +40°C | 32°F ... +104°F |
| Compatible with thermo washing machines and autoclavable < 135°C / 383°F (motor parts) | yes | yes |
| Noise level (motor only) acc. to ISO 11498 | 40 dBA | 40 dBA |

adeor® is a registered trademark of adeor Medical AG.

All rights of distribution, also by photocopy, reprint or storing and recovery in any kind of data processing device are reserved and require our written approval.

Copyright 2015 adeor Medical AG
01-0515-D/E
Printed in Germany
May 2015

Codman is a trademark of Johnson & Johnson, Inc.
adeor is not affiliated with any of the above.

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adeor reserves the right of
technical changes.

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adeor® ist eine eingetragene Marke der adeor Medical AG.
*Codman ist eine Marke der Johnson & Johnson, Inc.