



*2nd generation
available now*

CataRhex 3[®] easy

MAKING THE DIFFERENCE
WITH THE NEW PHACO POWER

ALL ADVANTAGES AT A GLANCE

NEW FEATURE

Larger control panel
Better visibility and legibility

Portable size
High mobility and space-saving

NEW FEATURE

Recessed grips
Easy moving and handling

NEW FEATURE

All plugs on front
Faster set-up of instruments during set-up

NEW FEATURE

Pneumatic push-pull principle
High speed cutting with pneumatic push-pull principle invented by Oertli

DirectAccess™
One button per function for simple operation

Precise flow control with peristaltic pump

easyPhaco® technology for safe and efficient emulsification

NEW FEATURE

Continuous Flow Cutter

- Continuously open port provides less traction force¹
- Constant cutting performance with up to 2'400 cpm thanks to 100 % quality checks²
- Precise single cut e.g. for iridectomy procedures

NEW FEATURE

Integrated air compressor
Independant application of cutter



¹ Compared to previous generation with SPS cutter
² 100 % final inspection with cutting test

Not available for sales in the US